



UA-704 Mini Manual Inflate



Digital blood pressure technology that fits in the palm of your hand!

The UA-704 Mini Manual Inflate is an innovative monitor that incorporates digital blood pressure technology and advanced features in a small, compact design. Included with this hand-held monitor is A&D's Irregular Heartbeat feature which allows those who have an irregular heartbeat a reliable method of taking their blood pressure at home.

A manual inflate monitor such as the UA-704 offers an economical solution to those who need an accurate way to monitor their blood pressure, but without the expense of an automatic inflate unit. In addition, this energy efficient monitor can perform approximately 2,000 measurements with just one 'AA' battery (included).

- Space-saver design ideal for travel
- Irregular Heartbeat feature
- Requires only one 'AA' battery (included)
- Available in medium or large cuff sizes (cuffs can also be purchased separately)
- Large-print instruction guide with illustrations and quick start card in English, French and Spanish
- Four week blood pressure logbook included

For questions about **LifeSource** monitors, call our toll free **LifeSource** HealthLine at 1-800-726-3364 and a specially trained representative will assist you, or visit us at www.LifeSourceOnline.com.

- New "Mini" size
- One button operation
- Professional accuracy
- Compact retail package
- Affordable and reliable
- Latex free



UA-704 Mini Manual Inflate

To obtain an accurate reading, it is important to use the correct cuff size. A cuff that is too small will result in a higher than expected reading and a cuff that is too large will result in a lower than expected reading. The UA-704 Mini Manual Inflate unit comes in two different cuff sizes to accommodate your customer's needs. This monitor is not compatible with the *LifeSource* small cuff.

	PART NUMBER	UPC NUMBER	CUFF SIZE
UNITED STATES AND CANADA	UA-704V	0-93764-60106-4	Medium 7.5 - 12.2" (19 - 31 cm)
	UA-704VL	0-93764-60107-1	Large 11.8 - 17.7" (30 - 45 cm)

Each unit comes with a trilingual instruction manual, quick start card, four week blood pressure logbook and warranty card.

Accessories

PART #	DESCRIPTION	PART #	DESCRIPTION
11-001	Latex Inflator Bulb	UA-280	Medium Blood Pressure Monitor Cuff: 7.5" - 12.2" (19-31cm)
		UA-281	Large Blood Pressure Monitor Cuff: 11.8" - 17.7" (30-45cm)
		UA-380	Medium <i>EasyCuff</i> ™ Feature: 7.5 - 14.2" (19 - 36 cm)

Technical Specifications

Dimensions	3.9" (l) x 1.2" (w) x 2.0" (h) (100mm x 31mm x 51mm)
Weight	2.6 oz. (73g) without batteries
Display	Digital: 10mm character height. Pressure and pulse displayed alternately
Pressurization	Manual, by inflator bulb
Deflation	Hand-held push button
Type	Oscillometric
Accuracy	Pressure: ±3mmHg or 2%, whichever is greater/Pulse: ±5%
Pressure Range	20mmHg to 280mmHg/Pulse: 40 to 200 pulses/minute
Power Source	One type AA alkaline batteries (1.5 volt). Battery included.
Battery Life	Approx. 2,000 measurements
Operating Environ.	50°F to 104°F (10°C to 40°C)
Storage Environ.	14°F to 140°F (-10°C to 60°C)

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9000 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Selling Unit (Box)	L x W x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
All UA-704 Models	3 ^{15/16} x 3 ^{15/16} x 7 ^{1/2}	0.93	0.07
Unit Case (10 units)	L x W x H (inches)	WT (lb.)	Case Cube (cubic feet)
All UA-704 Models	20 ^{1/2} x 8 ^{1/4} x 7 ^{7/8}	10.37	0.76

Ask your A&D Sales Representative about point-of-purchase floor and counter displays available for this blood pressure monitor unit.



1555 McCandless Drive
Milpitas, CA 95035
www.LifeSourceOnline.com

LifeSource™
Your source for a lifetime of health



LifeSource and EasyCuff are trademarks of A&D Medical. ©2003 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other *LifeSource* products.

- Blood Pressure Monitors
- Heart Rate Monitors
- Nebulizers
- Digital Thermometers
- Personal Scales
- Stethoscopes

For questions about *LifeSource* monitors, call our toll-free *LifeSource* HealthLine at 1-800-726-3364 or visit our web site at www.LifeSourceOnline.com

LITLS704 11/03