

## TALKING AUTO-INFLATION Blood Pressure Monitor UA-767T

**Fast. Easy. Accurate.**



**One button, fully automatic blood pressure monitor ideal for the visually impaired.**

A&D Medical adds the convenience of an audio announcement of measurement results to the popular One Step Auto-Inflation monitor. The cuff is automatically inflated at the touch of a button. When the measurement is complete, blood pressure and pulse are clearly heard as they are displayed, offering independent monitoring to those who prefer audio results in addition to digitally displayed results. The unit shuts off automatically after four announcements of the measurement or may be stopped at any time by pressing the start button. A stereo jack for headphones enhances privacy and convenience (headphones not included).

- One button operation
- Simultaneous readout of systolic/diastolic pressure and pulse rate
- Large, easy to read 16mm digital display
- Available in medium or large cuff sizes (cuffs can also be purchased separately)
- Large-print instruction guide with illustrations and quick start card in English, French and Spanish
- Four week blood pressure logbook included
- Study with Peer Review published in *Blood Pressure Monitoring*

For questions about LifeSource products, call us toll-free at 1-800-726-3364 and a specially trained representative will assist you, or visit us at [www.LifeSourceOnline.com](http://www.LifeSourceOnline.com).

- Clear Voice Announcement of Reading
- Clinically validate for accuracy
- One Button, Fully Automatic Operation
- Simple Preset Switch for Comfort
- Latex free

# UA-767T Talking Auto-Inflation

To obtain an accurate reading, it is important to use the correct cuff size. A cuff that is too small will result in a higher than expected reading and a cuff that is too large will result in a lower than expected reading. The UA-767T Talking Auto-Inflation unit comes in three different cuff sizes to accommodate your customer's needs.

PART NUMBER	UPC NUMBER	CUFF SIZE
UA-767T	0-93764-22767-7	Medium 9.4" - 14.2" (24 - 36cm)
UA-767TL	0-93764-60015-2	Large 14.2" - 17.7" (36 - 45cm)

## Accessories

PART NUMBER	UPC NUMBER	Description
UA-279	0-93764-01154-2	Small Cuff: 6.3" - 9.4" (16 - 24cm)
UA-280	0-93764-04607-0	Medium Cuff: 9.4" - 14.2" (24-36cm)
UA-281	0-93764-01069-9	Large Cuff: 14.2" - 17.7" (36-45cm)
UA-700CH	NA	Cuff Hose

Each unit comes with a trilingual instruction manual, quick start card, four week blood pressure logbook and warranty card.

## Technical Specifications

<b>Dimensions</b>	4.2" (l) x 6.4" (w) x 2.8" (h) (106mm x 163mm x 70mm)
<b>Weight</b>	10.9 oz. (308g) without batteries
<b>Display</b>	Digital: 16mm character height. Pressure and pulse displayed simultaneously
<b>Pressurization</b>	Automatic, using micropump
<b>Deflation</b>	Automatic exhaust
<b>Type</b>	Oscillometric
<b>Accuracy</b>	Pressure: ±3mmHg or 2%, whichever is greater/Pulse: ±5%
<b>Pressure Range</b>	20mmHg to 280mmHg/Pulse: 40 to 200 pulses/minute
<b>Power Source</b>	Four type AA alkaline batteries (1.5 volt). Batteries not included.
<b>Battery Life</b>	Approx. four months with one daily measurement
<b>Operating Environ.</b>	50°F to 104°F (10°C to 40°C)
<b>Storage Environ.</b>	14°F to 140°F (-10°C to 60°C)

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9000 and ISO 13485 registered company. Engineered in Japan. Made in China.

## Packaging and Shipping Specifications

<b>Selling Unit (Box)</b>	<b>L x W x H (inches)</b>	<b>WT (lbs.)</b>	<b>Individual Item Cube (cubic feet)</b>
UA-767T	5 <sup>1</sup> / <sub>16</sub> x 3 <sup>1</sup> / <sub>2</sub> x 7 <sup>11</sup> / <sub>16</sub>	1.30	0.08
UA-767TL	5 <sup>1</sup> / <sub>16</sub> x 3 <sup>1</sup> / <sub>2</sub> x 7 <sup>11</sup> / <sub>16</sub>	1.45	0.08
<b>Unit Case (10 units)</b>	<b>L x W x H (inches)</b>	<b>WT (lbs.)</b>	<b>Case Cube (cubic feet)</b>
UA-767T	18 <sup>1</sup> / <sub>2</sub> x 11 x 8 <sup>1</sup> / <sub>2</sub>	14.15	1.00
UA-767TL	18 <sup>1</sup> / <sub>2</sub> x 11 x 8 <sup>1</sup> / <sub>2</sub>	15.75	1.00



1756 Automation Parkway  
San Jose, CA 95131  
www.LifeSourceOnline.com

LIFE SOURCE®



LifeSource is a trademark of A&D Medical. ©2007 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Pedometers

For questions about LifeSource products, call us toll-free at 1-800-726-3364 or visit our web site at [www.LifeSourceOnline.com](http://www.LifeSourceOnline.com)

LIF:LS767T 12/07