


## Digital Thermometers



LifeSource digital thermometers provide accurate temperature readings at the touch of a button.

LifeSource digital thermometers provide a safe, convenient alternative to traditional glass mercury thermometers. A&D Medical offers three digital thermometers to meet different needs. Each thermometer is encased in durable plastic and are safe for all ages. A&D Medical's LifeSource thermometers provide readings in either 10 or 60 seconds.

The DT-703 Digital Thermometer provides 60 second measurements at an economical price point. The DT-704 Thermometer provides 10 second measurements with an extra large display and the DT-705 Thermometer provides 10 second measurements with a soft, flexible tip.

These thermometers have features that include a digital temperature readout, last measurement memory, audible beep when a measurement is complete and automatic shut off.

- Includes protective carrying case for safe storage
- Bilingual instruction guide with illustrations in English, and Spanish

For questions about LifeSource products, call us toll-free at 1-800-726-3364 and a specially trained representative will assist you.

- Accurate readings in Fahrenheit and Centigrade
- Mercury free
- Beeps when peak temperature reached
- Automatic shut-off
- Automatic last reading memory

# LifeSource Digital Thermometers

## DT-703 Digital Thermometer

- 60 second measurement time
- Reads Farenheit and Celsius
- Carrying case included
- UPC: 0-93764-60184-2



## DT-704 Fast-Read Thermometer

- 10 second measurement time
- Large display for easy reading
- Reads Farenheit and Celsius
- Carrying case included
- UPC: 0-93764-60185-9



## DT-705 Flex-Tip Thermometer

- 10 second measurement time
- Flexible tip for added comfort
- Reads Farenheit and Celsius
- Carrying case included
- UPC: 0-93764-60186-6



## Technical Specifications

<b>Range</b>	90.0° - 110.0° F (32.0° - 44.0° C)
<b>Accuracy</b>	± 0.2°F between 96.0° and 104.2° F; ± 0.1°C between 34.0° and 42.0° C. At standard room temperature of 64° - 82.4° F (18.0° - 28° C)
<b>Display</b>	0.01° F (0.01° C)
<b>Battery</b>	One 1.55V button size battery (LR41 or equivalent)
<b>Power Consumption</b>	Less than or equal to 0.09 mW
<b>Battery Life</b>	Approximately 3000 measurements
<b>Dimensions</b>	DT-703: 5.2" x 0.65" x 0.37" DT-704: 5.08" x 0.89" x 0.45" DT-705: 4.94" x 0.90" x 0.50"
<b>Weight</b>	Approximately 13 grams including battery
<b>Fast Read</b>	DT-703: Provides readings in approximately 60 seconds DT-704: Provides readings in approximately 10 seconds DT-705: Provides readings in approximately 10 seconds
<b>Beeper</b>	Beeps after peak temperature is reached.

## Packaging and Shipping Specifications

Selling Unit (Box)	L x W x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
DT-703	5.20" x 0.65" x 0.37	0.03	0.01
DT-704	5.08" x 0.89" x 0.45	0.03	0.01
DT-705	4.94" x 0.90" x 0.50"	0.03	0.01

Unit Case (20 units)	L x W x H (inches)	WT (lb.)	Case Cube (cubic feet)
DT-703	7.80" x 3.39" x 7.20"	1.98	0.11
DT-704	9.37" x 6.30" x 7.20"	2.65	0.12
DT-705	9.37" x 6.30" x 7.20"	2.87	0.12

Ask your A&D Sales Representative about countertop point-of-purchase displays available for these digital thermometers.



1756 Automation Parkway  
San Jose, CA 95131  
www.LifeSourceOnline.com



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes

For questions about LifeSource products, call us toll-free at 1-800-726-3364 or visit our web site at [www.LifeSourceOnline.com](http://www.LifeSourceOnline.com)