

DIGITAL WRIST Blood Pressure Monitor UB-521

Fast. Easy. Accurate.



- 90 Memory Recall
- Pressure Rating Indicator™
- Irregular Heartbeat Detection
- Displays average readings
- One Touch Button
- Ideal for travel

A&D Medical has applied its state-of-the-art measurement technology to monitoring blood pressure at the wrist with the push of a single button. This wrist monitor is smaller, lighter and more compact than other wrist monitors and is ideal for travel. Unlike other monitors that measure blood pressure during deflation of the cuff, the Travel Size wrist monitor measures blood pressure during cuff inflation for a faster, more comfortable reading. The automatic memory keeps your last 90 blood pressure measurements for effective tracking.

The average reading feature automatically calculates the “average” of the total readings stored in memory. The UB-521 also has the Irregular Heartbeat feature. This feature provides blood pressure and pulse rate measurements even when an Irregular Heartbeat occurs and alerts the user. A rigid carrying case is included to protect the monitor when not in use.

- One button operation
- Simultaneous readout of systolic/diastolic pressure and pulse rate
- Bilingual Instruction Manual with Blood Pressure Log Sheets (English and Spanish)
- Rigid carrying case included
- Includes two “AAA” batteries

For questions about LifeSource products, call us toll free at 1-800-726-3364 and a specially trained representative will assist you, or visit us at www.LifeSourceOnline.com.

Digital Wrist UB-521 Blood Pressure Monitor

The Digital Wrist blood pressure monitor fits easily and comfortably to one's wrist to provide an accurate reading. This compact and portable unit gives fast measurements and provides convenience to those looking for an easy way to monitor their blood pressure away from home.

PART NUMBER	UPC NUMBER	CUFF SIZE
UB-521	0-93764-60213-9	5.3" - 8.5" (13.5 - 21.5 cm)

Each unit comes with a bilingual instruction manual with blood pressure log sheets, quick start card, two AAA batteries and carrying case.

Technical Specifications

Dimensions	1.2" (d) x 2.4" (w) x 2.6" (h) (30 x 60 x 65 mm)
Weight	3.2 oz. (89 g) without batteries
Display	Digital: 9/7 mm character height Pressure and pulse displayed simultaneously
Memory	90 readings
Pressurization	Automatic, using micropump
Deflation	Constant air release-valve system
Type	Oscillometric
Accuracy	Pressure: ± 3 mmHg or $\pm 2\%$, whichever is greater; Pulse: $\pm 5\%$
Pressure Range	20mmHg to 280 mmHg; Pulse: 40 to 200 pulses/minute
Power Source	Two type AAA alkaline batteries (1.5 volt). Batteries included.
Battery Life	Approx. 400 measurements
Operating Environ.	50° to 104°F (10° to 40°C)
Storage Environ.	14° to 140°F (-10° to 60°C)

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9001 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Individual Unit (Box)	W x D x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
UB-521 models	4.6 x 3.52 x 6.12	0.71	0.06
Unit Case (10 units)	W x D x H (inches)	WT (lb.)	Case Cube (cubic feet)
UB-521 models	18 x 10 x 17	8.4	0.73



1756 Automation Parkway
San Jose, CA 95131
www.LifeSourceOnline.com

LIFE SOURCE®



Pressure Rating Indicator is a trademark of A&D Medical. Patent Pending. LifeSource is a registered trademark of A&D Medical. ©2008 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Pedometers

For questions about LifeSource products, call us toll-free at 1-800-726-3364 or visit our web site at www.LifeSourceOnline.com